



Med-Cut Scissors

with micro serrated blades



Med Surgical Instruments



Med-Cut Scissors

Precision Performance Every Time

Precision Performance Means *Razor Sharp*

- Microtome technology
 - Precision ground for the sharpest cutting edge


Every Time Means *Consistent High Performance*

- Durotip® TC technology
 - Tungsten carbide blades
- Smoothing technology
 - Low-wear internal blade surface
- Low-reflection coating
 - Superior wear resistance over full length





Med-Cut Scissors




Microtome technology ensures ultimate cutting edge precision for both smooth and Med-cut blade options. The razor-sharp microtome blade edge offers unparalleled sharpness with minimal use of force.

	Delicate				Delicate	
	1/1	1/1	1/1	1/1	1/1	1/1
	MS210B	MS211B	MS252B	MS253B	MS283B	MS284B
	straight	curved	straight	curved	curved	curved
	4-3/8"	4-3/8"	6-3/4"	6-3/4"	6-1/2"	6-1/2"
Iris Dissecting Scissors			MAYO Dissecting Scissors		MAYO-LEXER Dissecting Scissors	



Med-Cut Scissors



MS257B
curved
4-1/2"



Strabismus
Dissecting Scissors



METZENBAUM
Dissecting Scissors



MS261B
curved
5-3/4"



MS263B
curved
7"



MS271B
curved
7"



MS265B
curved
8"



MS275B
curved
8"



MS277B
curved
9"

Delicate

Delicate

Delicate




Med-Cut Scissors

The particularly delicate shape permits smooth dissection even in very delicate tissue structures.

	<p>Delicate</p>  <p>MS273B curved 7"</p>		 <p>MS259B curved 5-3/4"</p>		 <p>MS267B curved 9"</p>	<p>Delicate</p>  <p>MS279B curved 10-1/4"</p>	<p>Delicate</p>  <p>MS281B curved 11"</p>
 <p>TOENNIS-ADSON Dissecting Scissors</p>		 <p>BABY-METZENBAUM Dissecting Scissors</p>		 <p>NELSON-METZENBAUM Dissecting Scissors</p>			



Med-Cut Scissors



Smooth Blade Option

Item No.	Description
MS210 B	Iris Dissecting Scissors, delicate, straight, 4-3/8"
MS211 B	Iris Dissecting Scissors, delicate, curved, 4-3/8"
MS252 B	MAYO Dissecting Scissors, straight, 6-3/4"
MS253 B	MAYO Dissecting Scissors, curved, 6-3/4"
MS257 B	Strabismus Dissecting Scissors, curved, 4-1/2"
MS259 B	BABY-METZENBAUM Dissecting Scissors, curved, 5-3/4"
MS261 B	METZENBAUM Dissecting Scissors, curved, 5-3/4"
MS263 B	METZENBAUM Dissecting Scissors, curved, 7"
MS265 B	METZENBAUM Dissecting Scissors, curved, 8"
MS267 B	NELSON-METZENBAUM Dissecting Scissors, curved, 9"
MS271 B	METZENBAUM Dissecting Scissors, delicate, curved, 7"
MS273 B	TOENNIS-ADSON Dissecting Scissors, delicate, curved, 7"
MS275 B	METZENBAUM Dissecting Scissors, delicate, curved, 8"
MS277 B	METZENBAUM Dissecting Scissors, delicate, curved, 9"
MS279 B	NELSON-METZENBAUM Dissecting Scissors, delicate, curved, 10-1/4"
MS281 B	NELSON-METZENBAUM Dissecting Scissors, delicate, curved, 11"
MS283 B	MAYO-LEXER Dissecting Scissors, delicate, curved, 6-1/2"
MS284 B	MAYO-LEXER Dissecting Scissors, curved, 6-1/2"



Med-Cut Scissors



Med-cut micro serration technology prevents tissue from slipping out during the cutting process.

Delicate



MS211WB
curved
4-3/8"



Iris
Dissecting Scissors,
with Micro Serrated Blades



MS253WB
curved
6-3/4"



MAYO
Dissecting Scissors,
with Micro Serrated Blades

Delicate



MS283WB
curved
6-1/2"




MAYO-LEXER
Dissecting Scissors,
with Micro Serrated Blades



MS284WB
curved
6-1/2"



Med-Cut Scissors



1/1

MS257WB
curved
4-1/2"



Strabismus
Dissecting Scissors,
with Micro Serrated Blades



1/1

MS261WB
curved
5-3/4"



METZENBAUM
Dissecting Scissors,
with Micro Serrated Blades



1/1

MS263WB
curved
7"



1/1

MS271WB
curved
7"



1/1

MS265WB
curved
8"



1/1

MS275WB
curved
8"



1/1

MS277WB
curved
9"



1/1

Delicate




Med-Cut Scissors

The particularly delicate shape permits smooth dissection even in very delicate tissue structures.

	Delicate			Delicate	Delicate	Delicate	
	1/1		1/1	1/1	1/1	1/1	
	MS273WB curved 7"		MS259WB curved 5-3/4"		MS267WB curved 9"	MS279WB curved 10-1/4"	MS281WB curved 11"
TOENNIS-ADSON Dissecting Scissors, with Micro Serrated Blades		BABY-METZENBAUM Dissecting Scissors, with Micro Serrated Blades		NELSON-METZENBAUM Dissecting Scissors, with Micro Serrated Blades			



Med-Cut Scissors



Med-Cut Blade Option

Item No.	Description
MS211 WB	Iris Dissecting Scissors, delicate, with Micro Serrated Blades, curved, 4-3/8"
MS253 WB	MAYO Dissecting Scissors, with Micro Serrated Blades, curved, 6-3/4"
MS257 WB	Strabismus Dissecting Scissors, with Micro Serrated Blades, curved, 4-1/2"
MS259 WB	BABY-METZENBAUM Dissecting Scissors, with Micro Serrated Blades, curved, 5-3/4"
MS261 WB	METZENBAUM Dissecting Scissors, with Micro Serrated Blades, curved, 5-3/4"
MS263 WB	METZENBAUM Dissecting Scissors, with Micro Serrated Blades, curved, 7"
MS265 WB	METZENBAUM Dissecting Scissors, with Micro Serrated Blades, curved, 8"
MS267 WB	NELSON-METZENBAUM Dissecting Scissors, with Micro Serrated Blades, curved, 9"
MS271 WB	METZENBAUM Dissecting Scissors, delicate, with Micro Serrated Blades, curved, 7"
MS273 WB	TOENNIS-ADSON Dissecting Scissors, delicate, with Micro Serrated Blades, curved, 7"
MS275 WB	METZENBAUM Dissecting Scissors, delicate, with Micro Serrated Blades, curved, 9"
MS277 WB	METZENBAUM Dissecting Scissors, delicate, with Micro Serrated Blades, curved, 9"
MS279 WB	NELSON-METZENBAUM Dissecting Scissors, delicate, with Micro Serrated Blades, curved, 10-1/4"
MS281 WB	NELSON-METZENBAUM Dissecting Scissors, delicate, with Micro Serrated Blades, curved, 11"
MS283 WB	MAYO-LEXER Dissecting Scissors, delicate, with Micro Serrated Blades, curved, 6-1/2"
MS284 WB	MAYO-LEXER Dissecting Scissors, with Micro Serrated Blades, curved, 6-1/2"



Med-Cut Scissors



All rights reserved. Technical alterations are possible. This leaflet may be used for no other purposes than offering, buying and selling of our products. No part may be copied or reproduced in any form. In the case of misuse we retain the rights to recall our catalogs and price list and to take legal action.

Copyright 2012-14 Med Surgical. All rights reserved.



MED-SURGICAL INSTRUMENTS

724, Thornhill Drive
Coliegeville, PA 19426.
TEL: +1 610-983-3073

www.medsurgicalusa.com